


# PURINOZ

## Nasal Aspirator (Bean Style) BN-303N

1. The users can safely suck the mucus out from baby's nose through tube due to they can adjust the strength of suction manually by themselves.
2. Although baby's moving during aspiration, the nasal aspirator won't hurt baby's nasal mucosa due to circular-designed nose tip made of soft silicone.
3. Due to petal-type mouth piece, it can prevent the mouth piece from going too deep into user's mouth.
4. Due to the tube, mouth piece and nose piece are made of transparent silicone, the user can feel at ease during aspiration.
5. Due to unique design of anti-backflow device, the nasal aspirator can prevent the user from sucking the mucus through tube during aspiration.
6. It's easy to use and hold by user due to it's unique design of bean-shaped collection cup.



+ Bean shape + Easy to hold  
+ Easy clean + Easy to use  
+ Exchangeable Tips  
Material: Silicone, PP  
MADE IN TAIWAN

