

PURINOZ

Nasal Aspirator(Cricket Style) BN-506

The soft silicone nose cleaner can ease user's worry and protect baby's nasal mucosa



1. Nozzle and bulb syringe are made from non-toxic silicone and can withstand high temperatures.
2. Nozzle is made of soft silicone to avoid hurting baby's nasal mucosa.
3. Allow user to hold tightly due to its moderate size.
4. Press the bulb syringe from the left and right at the same time.



0+
mos

PAHs

- + Anti-backflow
- + Special dual way valve design
- + Food-grade silicone
- + High temperature resistance
- + Applicable in both nose and ear

Material:

Nozzle: Silicone/Bulb Syringe: Silicone/Tube: Polycarbonate

MADE IN TAIWAN



KenHealth®
www.kenhealth.biz