

SAFE BABY BATH THERMOMETER

BT08H

Baby Seahorse Bath Thermometer

Features

1. LED Screen.
2. Easy to Read.
3. Waterproof.
4. Hot Temperature Alarm (39°C).
5. Clock with Timing Display Function.
6. Changeable Battery.



0+
mos.

FDA

CE

ROHS

Specification

Measuring Range: 0° C to 50° C / 32° F to 122° F

Resolution: 1° C / 1° F

Accuracy: +/-1C during 25 to 50 degree C, Other Range is +/-2 Degree C

Module Dia: (Dia) 65 x (H) 50mm

Size: 138 x 90 x 51mm

C/F can be choice before produce.

Battery: 1 cell battery (1pc LR2032)

Battery Life: About one year (use one hour everyday)

Waterproof.

Thermometer Alarm Value (Default: 39.0° C).

Auto Off: Automatically into sleep mode after one hour

Material: Seahorse - PVC / Thermometer Module: ABS +PMMA

