

# Digital Thermometer

## TM-08K Digital Thermometer (Animal Head)

### Features

1. Idea for use by every family member.
2. For oral, rectal or underarm use.
3. Memory recall last temperature reading.
4. Automatic power off.
5. Fever alarm.
6. Beeper signal when temperature is ready to read.
7. Guaranteed clinically accurate.



0+  
mos.

FDA

CE

ROHS

### Specification

10 seconds reading.

Measuring Range: 32.0° C ~ 42.9° C (89.6° C ~ 109.2° F)

Accuracy: +/-0.1° C (+/-0.2° F)

Display: 3 1/2 LCD Digits (Unit of Display : 0.1° C / 0.1° F)

Battery: LR41 Alkaline

Battery Life: About 250 hours or 1000 times.

Weight: Appro. 16g

Unit Size: 151mm (L) x 24.5mm (W) x 18.5mm (H)

Waterproof.