

Digital Thermometer

TM-08K

Digital Thermometer (Economic)

Features

1. Idea for use by every family member.
2. For oral, rectal or underarm use.
3. Memory recall last temperature reading.
4. Automatic power off.
5. Fever alarm.
6. Beeper signal when temperature is ready to read.
7. Guaranteed clinically accurate.



Specification

With Soft Probe.

Measuring Range: 32.0° C ~ 42.9° C (89.6 ~ 109.2° F)

Measuring Time: About 10 seconds

Accuracy: $\pm 0.1^{\circ}$ C ($\pm 0.2^{\circ}$ F) for 35.0 ~ 42.0° C (95.0 ~ 107.6° F);
 $\pm 0.2^{\circ}$ C ($\pm 0.4^{\circ}$ F) for remaining range at an ambient of 22° C (71.6° F).

Display: 3 1/2 LCD digits (Unit of display: 0.1° C / 0.1° F)

Battery: LR41 Alkaline

Battery Life: About 250 hours or 1000 times.

Weight: 12g

Unit Size: 130 x 22 x 11mm

Waterproof.

