

Digital Thermometer

TM-09

Digital Thermometer (Flexible Tip)

Features

1. Idea for use by every family member.
2. For oral, rectal or underarm use.
3. Memory recall last temperature reading.
4. Automatic power off.
5. Fever alarm.
6. Beeper signal when temperature is ready to read.
7. Guaranteed clinically accurate.



Specification

Measuring Range: 32.0° C ~ 42.9° C (89.6° C ~ 109.2° F)

Accuracy: +/-0.1° C (+/-0.2° F)

60 Seconds reading.

Display: 3 1/2 LCD digits (Unit of display: 0.1° C / 0.1° F)

Battery: LR41 alkaline

Battery Life: About 250 hours or 1000 times.

Weight: Appro. 16g

Unit Size: 151mm (L) x 24.5mm (W) x 18.5mm (H)

Waterproof.

