

Digital Thermometer

TM-12F Digital Thermometer (Large Screen)

Features

1. Idea for use by every family member.
2. For oral, rectal or underarm use.
3. Memory recall last temperature reading.
4. Automatic power off.
5. Fever alarm.
6. Beeper signal when temperature is ready to read.
7. Guaranteed clinically accurate.



Specification

With Soft Probe.

C & F Selectable Before Production.

Measuring Range: 32.0° C ~ 42.9° C (89.6 ~ 109.2° F)

Measuring Time: 8 second

Accuracy: +/-0.1° C

Battery: 1 x CR1620

Battery Life: About 250 hours

Fever alarm.

Memory for last measurement.

Waterproof.

Automatic power off.

Low battery indication.