

Digital Thermometer

TM-21 Digital Thermometer (Red & Green Backlight)

Features

1. Idea for use by every family member.
2. For oral, rectal or underarm use.
3. Memory recall last temperature reading.
4. Automatic power off.
5. Fever alarm.
6. Beeper signal when temperature is ready to read.
7. Guaranteed clinically accurate.



Specification

Fever Glow.

Red / Green backlight after measurement.

C & F Selectable Before Production.

Measuring Range: 32.0° C ~ 42.9° C (89.6 ~ 109.2° F)

Measuring Time: 8 second

Accuracy: +/-0.1° C

Battery: 1 x CR2032

Battery Life: About 250 Hours

Fever Alarm.

Memory for last measurement.

Waterproof.

Automatic power off.

Low battery indication.

