

# PURINOZ



Electronic Nasal Aspirator 50KPA EN-03

Washable Nasal Aspirator- Gently clears  
stuffy nose in seconds



## Features :

- Power yet gentle suction
- Detachable collection cup set
- Fully immersible for rinse cleaning
- Cup set and nozzle able for boil-disinfection
- Unique design to wash inner pump
- 3 color coded soft silicone nozzles
- Zipper pouch included
- 2 X AA batteries required



## Specifications:

- Power: 2 X AA batteries ( not included)
- Suction: Approx. 50kPa.
- Size: 17 X 4.2 X 4.7 cm.
- Weight: 145g (exclude battery)
- Made in Taiwan



**KenHealth**<sup>®</sup>  
www.kenhealth.biz