

# Electronic Thermometer

LCT-200

Ear / Ambient thermometer

## Features

1. 2 in 1 Ear / Ambient thermometer.
2. One-second ear measurement.
3. Automatically stores 10 values.
4. Larged LCD screen for easy reading.
5. Over 1 minute idles, the power will turn off automatically for power saving.
6. Alert when the battery is low.
7. Ergonomically-shaped design.
8. Backlighted.
9. Beep sound when complete.
10. Fever alarm beeper, the light will turn green into orange color.
11. Easy to carry out with storage box.



Accuracy:  $\pm 0.2^{\circ}$  C

Operating Ambient Temperature Range:

$10^{\circ}$  C -  $40^{\circ}$  C ( $50^{\circ}$  F ~  $104^{\circ}$  F), Humidity Between 0% ~ 95%.

Displayed Temperature Range:  $0^{\circ}$  C -  $50^{\circ}$  C ( $32^{\circ}$  F ~  $122^{\circ}$  F)

Weight: 55g (with battery)

Memory: Auto-memorized the last 10 measured temperatures

Battery Type: CR2032 Lithium Cell (\*1)

Storage Box.

MADE IN TAIWAN

