

Electronic Thermometer

LCT-600

Ear / Ambient thermometer

Features

Features

1. 3 in 1 Ear / Ambient thermometer.
2. One-second ear measurement.
3. Automatically stores 10 values.
4. Larged LCD screen for easy reading.
5. Over 1 minute idles, the power will turn off automatically for power saving.
6. Alert when the battery is low.
7. Ergonomically-shaped design.
8. Backlighted.
9. Beep sound when complete.
10. Fever alarm beeper, the light will turn green into orange color.
11. Easy to carry out with storage box.



Accuracy: $\pm 0.2^{\circ}$ C

Operating Ambient Temperature Range:

10° C - 40° C (50° F ~ 104° F), Humidity Between 0% ~ 95%.

Displayed Temperature Range: 0° C - 50° C (32° F ~ 122° F)

Weight: 63g (with battery)

Memory: Auto-memorized the last 10 measured temperatures

Battery Type: CR2032 Lithium Cell (*1)

Storage Box.

MADE IN TAIWAN

