

# Fever Patch Sticker

TN-14

## Features

1. Efficiency, Continuity, Safety, Accuracy, and Hygiene.
2. Effectively indicates baby's temperature in both Fahrenheit & Centigrade.
3. Colored display.
4. Temperature measurement goes up by  $0.5^{\circ}$  C from  $35^{\circ}$  C to  $40^{\circ}$  C.
5. Easy to read.
6. Comfortable and simple to use.
7. Re-useable.
8. 12 patterns for choice.



0+  
mos.

## Specification

Temperature Scale :

35 - 35.5 - 36 - 36.5 - 37 - 37.5 - 38 - 38.5 - 39 - 39.5 -  $40^{\circ}$  C

Size of Product : 9.5cm x 2.5cm

Materials : Liquid Crystal Ink

