

Forehead Thermometer Strip (2 in 1 C&F) TN-07



Features

1. Efficiency, Continuity, Safety, Accuracy, and Hygiene.
2. Effectively indicates baby's temperature in both Fahrenheit & Centigrade.
3. Colored display.
4. Temperature measurement goes up by 1° C from 35° C to 40° C.
5. Easy to read.
6. Comfortable and simple to use.
7. Re-useable.
8. 12 patterns for choice.



0+
mos.

Specification

Temperature Scale:

- (1) Top > 95 - 96.8 - 98.6 - 100.4 - 102.2 - 104° F /
(2) Bottom > 35 - 36 - 37 - 38 - 39 - 40° C

Size of Product : 9.5cm x 2.5cm

Materials : Liquid Crystal Ink

