

Forehead Thermometer Strip (Popular) TN-27N

Features

1. Efficiency, Continuity, Safety, Accuracy, and Hygiene.
2. Effectively indicates baby's temperature in both Fahrenheit & Centigrade.
3. Colored display.
4. Temperature measurement goes up by 0.5° C from 35° C to 40° C.
5. Easy to read.
6. Comfortable and simple to use.
7. Re-useable.
8. 12 patterns for choice.



0+
mos.

Specification

Temperature Scale :

35 - 35.5 - 36 - 36.5 - 37 - 37.5 - 38 - 38.5 - 39 - 39.5 - 40° C

Size of Product : 9.5cm x 2.5cm

Materials : Liquid Crystal Ink

