

BABYFEAST

Manual Breast Pump (2 in 1) PBA-2201

MATERIAL USED :

Pump unit - ABS / Suction cap and base - PP

Suction pump - Silicone / Anti-dust breast shield cover - PP

Massage pad - Silicone / Funnel - PP

Anti-flow valve - Silicone / Milk collecting bottle - PP

Silicone tube - Silicone / Adapter - ABS / USB cable - PU

3 way connector - PP



0+
mos.

BPA
FREE

ISO 9001
ISO 13485

MATERIAL

Anti-dust breast shield cover-PP / Massage pad-Silicone

Anti-flow valve-Silicone / Milk collecting bottle-PP

Pump unit-ABS / Adapter-ABS / USB cable-PU

Suction kit Funnel kit Bottle stand-pp

