

# BABYFEAST

## Dexterous Electric Breast Pump SM-610

### FEATURES

Special designed for breastfeeding mothers, which is the best tool to express and store your milk.

One-hand operation, easy and convenient for a mother to express and relieve breast engorgement successfully Compact & light weight exquisite pump, convenient for outdoor activities



ISO 9001  
ISO 13485

### MATERIAL USED :

- Suction unit and battery compartment-ABS
- Silicone pad-silicone / Funnel-pp / Gasket-silicone
- Pump cleaning connector-pp / Bottle-pp
- Teat-silicone / Cap-pp

